



PRODUCT CATALOG FOR LAB



Dental 
Save Money. Business Better!



Guangzhou YiLinK Medical Equipment Technology Co.,Ltd
Tel: 020-32208271
Web: WWW.yilink-medical.com
Fax: 020-32208069
Add: Building F 405, Software Science Park, No.11 Caipin Road, Luogang Science City, Huangpu District, Guangzhou.



YiLiNK



PEEK

- Metal-free
- Low weight
- Extremely durable



reliable

YiLiNK



Material Characteristic:

Thermal Expansivity (25-500 c)[10/K]	10.5±0.5
Bending Strength [MPa]	400±60
Chemical Solubility [ug/cm ²]	≤ 40
Density[g/cm ³]	2.47±0.05

Highly compatible with excellent performance

PEEK is suitable for the digital manufacture of fixed and removable dentures, such as implant-supported superstructures, attachment or screw-retained restorations and also crowns and bridges including secondary and telescope crowns.

Recommended Indications

Fixed and removable dentures

Colors

Natural White Yellow Pink



Size


98*(10/12/14/16/18/20/22/25)mm

95*(10/12/14/16/18/20/22/25)mm

AG*(12/14/16/18/20/25)mm

Advantage

Density is near to natural teeth;
High resistance to wear, abrasion;
Easy to clean;



Hardness	95(Shore"A" Scale)
Density	0.92 g/cm³
Ignition point	298°C
melting point	114°C
viscosity	1790 CS
Volume shrinkage (natural point to room temperature)	6.24% typical
Ash content (lost wax casting)	0.0082%
Bending modulus	10.09
Coefficient of thermal expansion	1.09x10 ⁻⁴



Indications

Wax model

Calibration


Color



Blue

Red

White



wax

Highlights

- 1.Easier to process – Low resin content,cutting without bonding milling bur.
- 2.More environmentally friendly – No harmful dust and odour.Better match .
- 3.Better effect – No deformation after cutting, no burrs.
- 4.Higher success rate – Small thermal expansivity, greatly improve the success rate.