

Q L O N E  
D E N T A L P R O

# Qlone Dental Pro

A 'New' iPhone App  
for 3D Face Scanning



For information

0536 778 90 30  
info@onexdental.com

WHAT IS

# QLONE DENTAL PRO



- New smartphone-based **3D facial scanner** for dentistry
- From the creators of the popular Qlone app
- Designed specifically for **dental applications**

KEY FEATURES OF

# QLONE DENTAL PRO



- Up to **16x improvement** in **mesh quality** over original Qlone app
- **Detailed scans** in the teeth area for **alignment with intraoral scans**
- **Streamlined UI** for quick face scanning



# SOFTWARE COMPATIBILITY

Works with major dental CAD software:



BlueSkyBio

3shape 

exocad

MEDIT

Supports **STL**, **PLY**, and **OBJ** file formats



# HOW DOES THE APP WORK?



The scanning process is as follows:

- 1 Uses smartphone camera and AI algorithms
- 2 Guides users through capture process
- 3 Automatically optimizes resulting 3D model
- 4 Faster processing from scan to export

# ACCESSIBILITY & PRICING:



Available for **iOS devices**



**\$200 USD** per year



**Unlimited scans** and **exports**  
with **premium access**



**No automatic** renewal



# PRIVACY & SUPPORT:

- **Strict privacy policy** - no patient data stored on servers
- "Ask An Expert" feature for **software integration support**
- **Regular updates** promised, including **direct integrations** with dental software



# ADDITIONAL FEATURES

The Qlone app has other capabilities too:


- ✓ **Import** and **optimize** your own **3D models**
- ✓ **Automatically cropped** around **face area**
- ✓ **No need** for **manual cleaning** of 3D models
- ✓ **Potential replacement** for retired Bellus3D app



# WHAT IS THE POTENTIAL IMPACT OF THIS APP ON THE DENTAL PRACTICE?

- ✓ Enhanced **diagnosis capabilities**
- ✓ Improved **treatment visualization**
- ✓ Boosted **patient communication**
- ✓ More **accessible 3D facial scanning** for a broader range of practices





Will you be  
downloading the  
Qlone Dental  
Pro app?



For information

0536 778 90 30  
info@onexdental.com