

DIRECT 3D PRINT RETAINERS

Precision, Durability, and Comfort in Every Print

SprintRay Retainer resin is engineered to deliver high-quality, custom-fit retainers with exceptional durability and clarity. Designed exclusively for use with SprintRay Pro S and Pro 2 3D printers, this resin ensures that your patients receive retainers that not only fit perfectly but also stand the test of time.

This resin is Class I MDR certified and part of a fully validated biocompatible direct print retainer system, which includes free AI Retainer Design. With this breakthrough, dental professionals can look forward to producing high-quality retainers with precision and ease.



SprintRay Retainer
500g
SRI-0202161

FLEXURAL STRENGTH

RETAINER RESIN > 40 MPa*

FLEXURAL MODULUS

RETAINER RESIN > 1000 MPa*

TENSILE STRENGTH

RETAINER RESIN > 1000 MPa**

ULTIMATE TENSILE STRENGTH

RETAINER RESIN > 35 MPa**



KEY FEATURES

- ✓ **High Strength and Durability**
Long-lasting, durable retainers for daily protection.
- ✓ **Superior Fit and Comfort**
Precisely printed for a comfortable fit.
- ✓ **Crystal Clear Aesthetics**
Nearly invisible retainers with superior clarity.
- ✓ **Fast and Efficient Workflow**
Quick retainer production with free AI Retainer Design included.

*Testing method ASTM D790. **Testing method ASTM D638

SprintRay Retainer resin is compatible with the SprintRay Pro S 3D Printer and optimized for the SprintRay Pro 2 3D Printer. It is formulated to be cured exclusively with the SprintRay NanoCure curing unit.

info.es@sprinray.com Get in touch with our sales team to see how 3D printing can enhance your dental practice. Available for sale in the US.

www.sprinray.com